

# Tube Assortment

ORDER NO	mL	SIZE	PLAIN	WITH GEL ACTIVATOR	WITH EDTA K <sub>3</sub> SOLUTION	WITH SODIUM CITRATE	WITH LITHIUM HEPARINE	HEPARIN WITH EDTA K <sub>3</sub>
VWRT220	5 mL	12 x 25	○					
VWRT220	7 mL	12 x 100	○					
VWRT220	10 mL	16 x 100	○					
VWGA20	5 mL	12 x 25		○				
VWGA20	7 mL	12 x 100		○				
VWGA20	10 mL	16 x 100		○				
VWEX20	2 mL	12 x 75			○			
VWEX25	2,5 mL	12 x 75			○			
VWEX30	3 mL	12 x 75			○			
VWEX35	3,5 mL	12 x 75			○			
VWEX40	4 mL	12 x 75			○			
VWEX45	4,5 mL	12 x 75			○			
VWEX70	7 mL	12 x 100			○			
VWEX100	10 mL	16 x 100			○			
VWRT300	10 mL	12 x 75				○		
VWRT350	5,5 mL	12 x 75				○		
VWLR20	20 mL	12 x 75					○	
VWLR40	4 mL	12 x 75					○	
VWLR50	5 mL	12 x 100					○	
VWLR80	8 mL	12 x 100					○	
VWLR70	7 mL	12 x 100					○	
VWLR90	9 mL	16 x 100					○	
VWLR90	9 mL	16 x 100					○	
VWEP001								○
HOLDER								